




# DENDRIDIAG® IW

ATP tests for pulp and paper industry



## MEASURE THE QUANTITY OF BACTERIA IN YOUR WATER CIRCUIT IN UNDER 2 MINUTES

### MANAGE WATER RECYCLING

Control microbiological growth associated with water recycling.

### ENSURE PRODUCTION CONTINUITY

Limit the formation of slime (biofilm) leading to paper breaks, stains, discoloration...

### ENVIRONMENTAL IMPACT REDUCED

Master the use of treatment products and limit highly polluted wastewater (filamentous bacteria, biocide, bio dispensant...).

## MASTER YOUR WATER CIRCUIT



#### Monitor

the microbial quality of your water



#### Identify

the critical points in no time at all



#### Take action

immediately in a crisis situation



#### Optimize

your cleaning and disinfection process


## WHAT IS DENDRIDIAG®?



A complete kit to measure the quantity of bacteria in industrial water

◀ [Watch the video](#)


- Cooling water systems,
- Industrial processes,
- Electroplating baths...




Get the result in a second thanks to the **web app!**

# A FAST AND RELIABLE ANALYSIS


## 1 FILTER the sample




## 2 EXTRACT the ATP



## 3 QUANTIFY the ATP



## 4 VALIDATE the analysis



## INDIVIDUAL SUPPORT

### SET UP

- Remote presentation of the kit
- Identification of needs
- Product shipping and scheduling your training date

1 week

### INDIVIDUAL SUPPORT

- Result analysis assistance
- Validation of the sampling points
- Confirmation of the thresholds

1 month

### ROLLOUT

- On-site or remote training
- Identification of the sampling points
- Definition of analysis frequency

1 year

### QUALITY CONTROL


- Device verification
- Customer feedback
- Forward planning

## KIT CONTENT



### The reagents

- 60 measurement tubes LUMITUBE
- 1 EXTRACTANT dropper
- 1 STANDARD dropper



### The consumables

- 60 sampling containers
- 60 10 ml syringes
- 60 0.45 µm filters



### The measuring device

- 1 luminometer KIKKOMAN SMART
- 1 microtube rack
- 1 USB flash drive

## NEED INFORMATION? CONTACT US



### Head office

5, avenue de l'Europe - Hélioparc  
34 830 CLAPIERS (FRANCE)

### Grand Est agency

8, rue Alfred Kastler  
Maison de l'entreprise  
54 320 MAXÉVILLE (FRANCE)



+ 33 (0)9 67 39 35 20



[contact@gl-biocontrol.com](mailto:contact@gl-biocontrol.com)



Join us on **LinkedIn**

**GLBIOCONTROL**

RAPID MICROBIAL DIAGNOSTIC

[WWW.GL-BIOCONTROL.COM](http://WWW.GL-BIOCONTROL.COM)