



DENDRIDIAG® IW

ATP tests for industrial laundry

MEASURE THE QUANTITY OF BACTERIA
IN YOUR WATER CIRCUIT IN UNDER 2 MINUTES

R.A.B.C. APPROACH

Proactively address microbiological deviations that may result in regulatory non-compliance.

SAVE

Optimize the dosage of laundry products, washing temperatures, load ratios...

ENVIRONMENTAL IMPACT REDUCED

Decrease water consumption and limit the release of harmful agents into the environment (AOX, phosphorus, detergents...).

MASTER YOUR WATER CIRCUIT



Monitor

the microbial quality of your water



Identify

the critical points in no time at all



Take action

immediately in a crisis situation



Optimize

your cleaning and disinfection process


WHAT IS DENDRIDIAG®?



A complete kit to measure the quantity of bacteria in industrial water


◀ **Watch the video**

- Cooling water systems,
- Industrial processes,
- Electroplating baths...




Get the result in a second thanks to the **web app!**

A FAST AND RELIABLE ANALYSIS



INDIVIDUAL SUPPORT




KIT CONTENT



The reagents

- 60 measurement tubes LUMITUBE
- 1 EXTRACTANT dropper
- 1 STANDARD dropper



The consumables

- 60 sampling containers
- 60 10 ml syringes
- 60 0.45 µm filters



The measuring device

- 1 luminometer KIKKOMAN SMART
- 1 microtube rack
- 1 USB flash drive

NEED INFORMATION? CONTACT US



Head office

5, avenue de l'Europe - Hélioparc
34 830 CLAPIERS (FRANCE)

Grand Est agency

8, rue Alfred Kastler
Maison de l'entreprise
54 320 MAXÉVILLE (FRANCE)



+ 33 (0)9 67 39 35 20



contact@gl-biocontrol.com



Join us on **Linkedin**

GLBIOCONTROL

RAPID MICROBIAL DIAGNOSTIC

WWW.GL-BIOCONTROL.COM